Surgery Data Manager

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The main settings area is titled '設定' (Settings) and includes the following sections:

- 出力ルート (ネットワークHDD/NAS)** (Output Path): A text field shows '/Volumes/Stella_8TB'. A button 'フォルダを選択' (Select Folder) is next to it. A yellow warning box below states: 'NASがマウントされていない場合は事前にFinderで接続してください。`/Volumes/` 以下のパスを指定します。' (If NAS is not mounted, please connect it in Finder in advance. Specify the path below ` /Volumes/ `).
- DICOM 送信先 (StellaDICOM 等)** (DICOM Destination):
 - Calling AE Title (自局): A text field with 'SURGERY'.
 - Called AE Title (相手局): A text field with 'STELLADICOM'.
 - Host: A text field with 'stelladicom.local'.
 - Port: A text field with '104'.
 - Modality: A dropdown menu with 'OT' selected.
 - Transfer Syntax: A dropdown menu with 'Implicit VR Little Endian' selected.
 - A button 'C-ECHO で疎通確認' (Check connectivity with C-ECHO) is located below these fields.
- 種別ごとの保存先フォルダ** (Folder for each type):
 - A note: '各種別のデータは指定フォルダ配下に「患者フォルダ」を作って収納されます (例: 手術写真/2026-05-19_P0001_モモ_去勢術/)。絶対パス推奨。空欄や相対パスにすると「出力ルート」配下として解決されます。' (Data of each type is stored in the specified folder hierarchy by creating a 'Patient Folder' (e.g., Surgery Photo/2026-05-19_P0001_Momo_Castration). Absolute path is recommended. If empty or relative path is used, it will be resolved under the 'Output Path').
 - 麻酔モニター記録** (Anesthesia Monitor Record): A text field shows '/Volumes/Stella_8TB/麻酔記録'. A button '選択...' (Select...) is next to it.
 - 手術写真 (DICOM送信対象)** (Surgery Photo (DICOM Transfer Target)): This section is partially visible at the bottom.

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